

## Product datasheet for **RR210284L3V**

### Ggt5 (NM\_019235) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ggt5 (NM_019235) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ggt5
Synonyms:	Ggta1; Ggtla1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_019235
ORF Size:	1716 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210284).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_019235.2</a> , <a href="#">NP_062108.2</a>
RefSeq Size:	2636 bp
RefSeq ORF:	1719 bp
Locus ID:	29566
UniProt ID:	<a href="#">Q5PQV0</a>
Cytogenetics:	20p12
Gene Summary:	involved in protecting the lungs from stress due to oxidation or other types of toxics [RGD, Feb 2006]



[View online »](#)