

## Product datasheet for RR210286L4

### Psmc1 (NM\_057123) Rat Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Psmc1 (NM_057123) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Psmc1
Synonyms:	s4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210286).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.



[View online >](#)



Locus ID: 117263  
UniProt ID: [P62193](#)  
Cytogenetics: 6q32  
Gene Summary: may play a role in ATP-dependent RNA/DNA unwinding [RGD, Feb 2006]