

Product datasheet for **RR210298L4V**

Casc3 (NM_147144) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Casc3 (NM_147144) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Casc3
Synonyms:	Btz; Mln51
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_147144
ORF Size:	2097 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210298).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_147144.1 , NP_671485.1
RefSeq Size:	2259 bp
RefSeq ORF:	2100 bp
Locus ID:	259170
UniProt ID:	Q8K3X0
Cytogenetics:	10q31
Gene Summary:	nucleocytoplasmic shuttling protein, component of the dendritic mRNA transport particles in hippocampal neurons [RGD, Feb 2006]



[View online »](#)