

## Product datasheet for **RR210354L4V**

### F2r11 (NM\_053897) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	F2r11 (NM_053897) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	F2r11
Synonyms:	Par-2; Par2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053897
ORF Size:	1191 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210354).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_053897.2</a> , <a href="#">NP_446349.1</a>
RefSeq Size:	3011 bp
RefSeq ORF:	1194 bp
Locus ID:	116677
UniProt ID:	<a href="#">Q63645</a>
Cytogenetics:	2q12
Gene Summary:	G-coupled protein receptor that is activated by coagulation factor II and other serine proteases; may play a role in regulation of exocrine secretion [RGD, Feb 2006]



[View online »](#)