

Product datasheet for **RR210365L3V**

Slc38a7 (NM_001003705) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc38a7 (NM_001003705) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Slc38a7
Synonyms:	MGC105833
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001003705
ORF Size:	1389 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210365).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001003705.2 , NP_001003705.1
RefSeq Size:	1840 bp
RefSeq ORF:	1392 bp
Locus ID:	291840
UniProt ID:	Q6JWR2
Cytogenetics:	19p13
Gene Summary:	Mediates sodium-dependent transport of amino acids, preferentially L-glutamine. [UniProtKB/Swiss-Prot Function]



[View online »](#)