

## Product datasheet for **RR210390L3V**

### Atad3a (NM\_001034922) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Atad3a (NM_001034922) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Atad3a
Synonyms:	Atad3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001034922
ORF Size:	1773 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210390).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001034922.1</a> , <a href="#">NP_001030094.1</a>
RefSeq Size:	2416 bp
RefSeq ORF:	1776 bp
Locus ID:	298682
UniProt ID:	<a href="#">Q3KRE0</a>
Cytogenetics:	5q36



[View online »](#)

**Gene Summary:**

Essential for mitochondrial network organization, mitochondrial metabolism and cell growth at organism and cellular level. May play an important role in mitochondrial protein synthesis. May also participate in mitochondrial DNA replication. May bind to mitochondrial DNA D-loops and contribute to nucleoid stability. Required for enhanced channeling of cholesterol for hormone-dependent steroidogenesis (By similarity).[UniProtKB/Swiss-Prot Function]