

## Product datasheet for **RR210420**

### Tsc22d3 (NM\_031345) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Tsc22d3 (NM\_031345) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Tsc22d3  
**Synonyms:** Dsipi; Gilz  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR210420 representing NM\_031345  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGAACACCGAAATGTATCAGACGCCCATGGAGGTGGCGGTCTATCAACTGCACAATTTCTCCATCTCCT  
TCTTCTCTTCTGCTTGGAGGGGATGTGGTTCCGTTAAACTGGATAACAGTGCCTCTGGAGCCAGCGT  
GGTGGCCCTAGACAACAAGATTGAGCAGGCCATGGATCTAGTGAAGAATCATTGATGTATGCTGTGAGA  
GAGGAGGTGGAGGTCCTAAAGGAGCAGATTTCGTGAGCTGGTTGAGAAGAACTCGCAGCTGGAGCGTGAGA  
ATACCCCTCTGAAGACGCTGGCGAGCCCAGAGCAGCTGGAGAAGTTCCAGTCCCAGCTGAGCCCCGAAGA  
GCCAGCTCCTGAAGCCCCAGAAACCCAGAAAGCCCTGGCGGTTCTGCGGTG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR210420 representing NM\_031345  
Red=Cloning site Green=Tags(s)

MNTEMYQTPMEVAVYQLHNFSISFFSLLGGDVSVKLDNSASGASVVALDNKIEQAMDLVKNHLMYAVR  
EEVEVLKEQIRELVEKNSQLERENTLLKTLASPEQLEKFQSRLSPEEPAPEAPETPEAPGGSAV

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI



[View online »](#)



<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_031345.1</a> , <a href="#">NP_112635.1</a>
<b>RefSeq Size:</b>	1969 bp
<b>RefSeq ORF:</b>	405 bp
<b>Locus ID:</b>	83514
<b>UniProt ID:</b>	<a href="#">Q9EQZ1</a>
<b>Cytogenetics:</b>	Xq32
<b>MW:</b>	14.8 kDa
<b>Gene Summary:</b>	may play a role in urinary Na(+) and K(+) excretion [RGD, Feb 2006]