

Product datasheet for **RR210428L3V**

Bbox1 (NM_022629) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Bbox1 (NM_022629) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Bbox1
Synonyms:	Bbh; Bbox
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022629
ORF Size:	1161 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210428).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022629.1 , NP_072151.1
RefSeq Size:	1164 bp
RefSeq ORF:	1164 bp
Locus ID:	64564
UniProt ID:	Q9QZU7
Cytogenetics:	3q34
Gene Summary:	Catalyzes the formation of L-carnitine from gamma-butyrobetaine.[UniProtKB/Swiss-Prot Function]



[View online »](#)