

Product datasheet for **RR210471**

Caly (NM_138915) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Caly (NM_138915) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Caly
Synonyms:	Calcyon; Drd1ip
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RR210471 representing NM_138915 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTGAAACTAGGCTGCAGCTTCTCAGGGAAGCCAGGCAAGGAAACAGGAGACCAGGATGGGGCTGCCA
TGGACAGTGTGCCTCTGATCAGCCCCCTGGATGTCAGCCAGCTTCAGCCATCATTCCCTGACCAGGTGGT
CATCAAAACACAGACAGAATATCAACTTACCTCTGCGGACCAGCCAAAGAAGTTCGCAGATTTGGAGGGT
CAGAGGCTGGCCTGCAGCCACCCAGAGGAAGGACGCAGACTGCCCACTGCAAGGATGATTGCCTTTGCTA
TGGCACTTCTGGGCTGTGTGCTGATCATGTACAAGGCCATTTGGTATGACCAGTTCACCTGCCAGATGG
CTTCTACTTCGGCACAAGATCTGCACCCCACTGACCCTGGAGATGTACTACACCGAGATGGACCCCGAA
CGCCACCGCAGCATCCTGGCAGCCATCGGGCCTACCCACTGAGCCGCAAGCACGGCACTGAGATGCCCG
CCATCTGGGGAACAGCTACCGGGCCGCAAGGAGGAGCATAAGGGCACCCAGCTGCGATGACTGT
GTCCACTGCAGCTGCAGCTGCAGCGGCAGAGGGCAACGAGCCGTGAGGAAAGCCTTTGGACATGAGAGAG
AAAGAAGATCCGCAGAAGGCGGAGGATGTGCCGTCCCAGTCCCCAAGAAAGAAGATCCGCAGAAGGCGG
AGGATGTGCCGTCCCAGTCCCCAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR210471 representing NM_138915
Red=Cloning site Green=Tags(s)

MVKLGCSFSGKPGKETGDQDGAAMDSVPLISPLDVSQ LQPSFPDQVVIKTQTEYQLT SADQPKKFADLEG
 QRLACSHPEEGRRLPTARMI AFAMALLGCVLIMYKAIWYDQFTCPDGFLLRHKICTPLTLEMYYTEMDPE
 RHR SILAAIGAYPLSRKHGTMPAIWGN SYRAGKEEHKGTTPAAMTVSTAAAAAAAE GNEPSGKPLDMRE
 KEDPQKAEDVPSQSPKKEDPQKAEDVPSQSPK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja1652_f05.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_138915

ORF Size: 678 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq Size: 1050 bp

RefSeq ORF: 681 bp

Locus ID: 192349


UniProt ID: [P58821](#)

Cytogenetics: 1q41

MW: 24.7 kDa

Gene Summary: dopamine receptor 1 (D1R) of D1-like receptors interacting protein [RGD, Feb 2006]

Product images:



Circular map for RR210471