

## Product datasheet for **RR210494**

### **Il7r (NM\_001106418) Rat Tagged ORF Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Il7r (NM_001106418) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Il7r
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR210494 representing NM_001106418 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGATGGCTCTGGGTAGAGCTTTCGCTATAGTTTTCTGCTTACTTCAAGCTGCTTCTGGAGAAAGTGGCA  
ATGCCCAGGATGGAGATCTAGAGGATGCGGAACCAGATGATCACTCCTTCTGGTCCACAGCCAGCTGGA  
AGTGGATGGAAATCAGCACTCACTGACGTGTGCTTTTAATGACCCAGACATCAAACCCTAATCTGGAA  
TTTCAAATATGTGGGCTCTTCTAGGCATAGATTGCCTAACTCTTAATAAGCTACGAGAGATGATTTTA  
TAAAGACATCAAATCTTACTGATTGGTAACAGCAGTGTATGTGTGAAGCTTGGAAAAATGGATGTAAT  
TTGCAAAATTTTGGACATAAGCACAATAGTTAAACCTGAGGCGCCTTCTAACCTGAAAGTAGTTTATCGA  
AAAGAAGCAAATGACTTTTTGGTAACATTTAATACATCTCACTCAACAAAGAAATATGTAACAGCATTAA  
AGCATGATGTGGCTACCGCCAGAAAGGGGTGAAAGTAACTGGACGCATGTATATTTATCCATACAAG  
AACAAATCCTACAGAGAAAATAACAACAAAGCAGTGTATGAAATCAAAGTCCGATCCATTCCAAT  
CATGAATACTTCAAAGGCTTCTGGAGTGTGGAGTCCAAGTCTACCTTTGAAACACCAGATCCAAGT  
ATCAAGGGGGATGGGATCCGGTTTTGCCAAGTATCATCCTCTGAGTTTGTCTCTATGGTTTTGTGGT  
CATTTTAGCCCATGTGCTATGGAAAAAAGGATTAACCTGTTGTATGGCCTAGTCTTCTGATCATAAG  
AAAACCTCTGGAACAATATGTAAGAAGCAAAAAAGAATTTGAATGTGAGTTTCAATCCTGAAAGTTTCC  
TGGACTGCCAGATTCATGAGGTGAATGGCATTCAAGCCAGGGATGAAGTGGAAAGCTTCTGCAAAATGA  
TCTTCTCCACGGCTGGGGAGCTGGAGAAGCAGGGACATAGAGCAACTGTACACGGTCAAACCTGGCCA  
TCTGAGATTTAGGCAGCACACCAGAAACGTTCCGAAGAGAATCACCCCTAAGATGCCTGGCTAGAAATC  
TAAGCACATGCAATACCCTGCATTCTTCTCTAGGTCCTGACTACAGAGAAGGTGACGGAAATAG  
GTCTCATGTGTATCAAGACTTGTGCCAAGCTCCAGAAACACAAATGGCACTGTCCCTCAACCATTTCT  
CTCCAGTCAGGAATCCTGATACCAAGTTTCTCAAGGACAGCCATCTACTTCTCAGTATTGAATCAAG  
AAGAAGCATATGCACTATGTCTAGCTTTTACCAAAACAAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


**Protein Sequence:** >RR210494 representing NM\_001106418  
 Red=Cloning site Green=Tags(s)

MMALGRAFAIVFCLLQAASGESGNAQDGDLEDAEPDDHSFWCHSQLEVDGNQHSLTCAFNDPDIKTTNLE  
 FQICGALLGIDCLTLNKLREMYFIKTSKFLIGNSSVCVKLGKMDVICKILDISTIVKPEAPSNLKVVYR  
 KEANDFLVTFNTSHSTKKYVTALKHDVAYRPERGESNWITHVYLFHTRTTILQRKLQPKAVYEIKVRSIPN  
 HEYFKGFWEWSPSSTFETPDSKYQGGWDPVLPISIILLSLFSMVLVLAHVLWKKRIKPVVWPSLPDHK  
 KTLEQLCKKPKKNLNVSNPESFLDCQIHEVNGIQARDEVESFLQNDLPPRPGELEKQGHRAVHGANWP  
 SEISGSTPETFRRESPLRCLARNLSTCNTPAFLSSRSPDYREGDGNRSHVYQDLLPSSRNTNGTVQPFP  
 LQSGILIPVSQGPPISTSSVLNQEEAYVTMSSFYQNK


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



<b>ACCN:</b>	NM_001106418
<b>ORF Size:</b>	1371 bp
<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_001106418.1</a> , <a href="#">NP_001099888.1</a>
<b>RefSeq Size:</b>	3124 bp
<b>RefSeq ORF:</b>	1374 bp
<b>Locus ID:</b>	294797
<b>Cytogenetics:</b>	2q16
<b>MW:</b>	51.5 kDa