



Product datasheet for RR210498

Cd28 (NM_013121) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cd28 (NM_013121) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Cd28
Synonyms:	CD28RNA
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:




* The last codon before the Stop codon of the ORF



[View online »](#)

Plasmid Map:



ACCN: NM_013121

ORF Size: 654 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_013121.1](#), [NP_037253.1](#)

RefSeq Size: 1284 bp

RefSeq ORF: 657 bp

Locus ID: 25660

UniProt ID: [P31042](#)

Cytogenetics: 9q32

MW: 25.2 kDa

Gene Summary: costimulatory molecule involved in induction of T cell proliferation and initiation of the Th2 response; may be important in maintenance of basal IgE production; binds to B7-1 or B7-2 on antigen presenting cells (APCs) [RGD, Feb 2006]