

Product datasheet for **RR210508**

Rps24 (NM_031112) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rps24 (NM_031112) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Rps24
Synonyms: L33
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR210508 representing NM_031112
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAATGACACGGTAACCATCCGGACCAGGAAGTTCATGACAAACCGTCTGCTTCAGAGGAACAGATGG
TCATTGATGTCTTCATCCTGGGAAGGCCACAGTACCAAAGACAGAAATTCGGGAAAAGCTGGCCAAGAT
GTACAAAACCACACCAGATGTCATCTTTGTATTTGGATTCAGAACCACCTTTGGTGGAGGCAAGACAAC
GGCTTTGGCATGATCTATGATTCTTTAGATTATGCAAAGAAGAATGAGCCTAAACACAGACTTGGCAGAC
ATGGCCTTTATGAGAAGAAAAGACCTCCCGCAAACAGCGAAAAGAACGAAAGAACAAGAATGAAGAAGGT
CAGGGGGACTGCAAAGGCCAATGTTGGTGTGGCAAAAAGCCAAAGGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR210508 representing NM_031112
Red=Cloning site Green=Tags(s)

MNDTVTIRTRKFMNRLLRKQMVIDVLHPGKATVPKTEIREKLAKMYKTPDVIFVFGFRTHFGGKKT
FGFMIYDSLIDYAKKNEPKHRLARHGLYEKKTSRKQRKERNMCKVVRGTAKANVGAGKKPKE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_031112.1 , NP_112374.1
RefSeq Size:	466 bp
RefSeq ORF:	402 bp
Locus ID:	81776
UniProt ID:	P62850
Cytogenetics:	16p16
MW:	15.4 kDa
Gene Summary:	ribosomal protein subunit [RGD, Feb 2006]