

Product datasheet for RR210563L4

Armcx3 (NM_001014273) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Armcx3 (NM_001014273) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Armcx3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210563).
Restriction Sites:	Sgfl-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_001014273
ORF Size:	1137 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001014273.1 , NP_001014295.1
RefSeq Size:	1992 bp
RefSeq ORF:	1140 bp
Locus ID:	367902
UniProt ID:	Q5XID7
Cytogenetics:	Xq32
Gene Summary:	Regulates mitochondrial aggregation and transport in axons in living neurons. May link mitochondria to the Trak2-kinesin motor complex via its interaction with Miro and Trak2. Mitochondrial distribution and dynamics is regulated through Armxc3 protein degradation, which is promoted by PCK and negatively regulated by Wnt1. Enhances the Sox10-mediated transactivation of the neuronal acetylcholine receptor subunit alpha-3 and beta-4 subunit gene promoters.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RR210563L4