

Product datasheet for **RR210612L4V**

Prpmp5 (NM_172065) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Prpmp5 (NM_172065) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Prpmp5
Synonyms:	Prb1; Prpg2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172065
ORF Size:	885 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210612).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172065.2 , NP_742062.1
RefSeq Size:	1010 bp
RefSeq ORF:	888 bp
Locus ID:	257651
UniProt ID:	P10165
Cytogenetics:	4q42
Gene Summary:	member of a family of proline rich proteins that are induced in salivary gland by isoproterenol treatment and by feeding tannins [RGD, Feb 2006]



[View online »](#)