

Product datasheet for **RR210618L3V**

Nudcd2 (NM_001009621) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nudcd2 (NM_001009621) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Nudcd2
Synonyms:	RGD1307203
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001009621
ORF Size:	471 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210618).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001009621.2 , NP_001009621.1
RefSeq Size:	1569 bp
RefSeq ORF:	474 bp
Locus ID:	287199
UniProt ID:	Q5M823
Cytogenetics:	10q12
Gene Summary:	May regulate the LIS1/dynein pathway by stabilizing LIS1 with Hsp90 chaperone. [UniProtKB/Swiss-Prot Function]



[View online »](#)