

Product datasheet for **RR210649**

Tmem182 (NM_001109305) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tmem182 (NM_001109305) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Tmem182
Synonyms: RGD1563696
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR210649 representing NM_001109305
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGAGACTAAATGTTGCTATTTTCTTTGGGCTCTCTTCGGAGCTTTGGGGTTTTGTTGTTTTGGTAG
CGTTTGGATCGGACTATTGGCTTCTTGCACCGAAAGTGGGAGGTGTTCCAGGTGAACAGAATATAGAGAA
CATCACTTTTACCATGAAGGCTTCTTCTGGAGGTGCTGGTTCAAAGGCGTTGTGGAGGAGAACAACCTCC
AACATCTGGAAGTTCTGGTACACCAACCAGCCGCATCCAAGAACTGCACACACGCTTACCTGTCCCGT
ACCCCTTCATGAGGGGTGAGCACAACCTCCACCTCCTATGACTCCGCAATTTTACCGTGGCTTCTGGGC
GGTCTGCTGCTCCTGGGGTCTGGTGCCTTGACTGCAAGCTTCTGATCATCTGCGCAGCCCCCTTT
TCCAGTCACTTACTTTATAAAGCTGGAGGAGGCTCATACATTGCCTCAGGTGTCCTGTTTTCCCTGGTGG
TGATACTGTATGTCATCTGGGTCCAGGCAGTGTCTGACATGAAAAGCTACAGAGCCTCGAAAATGAAAGA
CTGCTGGGAGTTCACGCCTTCTGTTCTGTACGGCTGGTCGTTTTTCTGGCCCCGCTGGAGATTTTTTC
TCTCTGCTCGCTGGGCTACTTTTCTTGTGTTGGGTGGCATATCCAGACGCATCAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR210649 representing NM_001109305
 Red=Cloning site Green=Tags(s)

MRLNVAIFFGALFGALGVLLFLVAFGSDYWLLATEVGRCSGEQNIENITFHHEGFFWRCWFKGVVEENNS
 NIWKFWYTNQPPSKNCTHAYLSPYPMRGEHNSTSYDSAIIYRGFWAVLLLLLVVAALTAFLIICAAPF
 SSHLLYKAGGGSYIASGVLFSLVVILYIWWQAVSDMESYRASKMKDCWEFTPSVLYGWSFFLAPAGVFF
 SLLAGLLFLVVGWHIQTHH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI


Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

Plasmid Map:



ACCN: NM_001109305

ORF Size: 687 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001109305.1 , NP_001102775.1
RefSeq Size:	2708 bp
RefSeq ORF:	690 bp
Locus ID:	501129
Cytogenetics:	9q22
MW:	25.8 kDa