

Product datasheet for **RR210650**

Npas1 (NM_001107479) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Npas1 (NM_001107479) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Npas1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RR210650 representing NM_001107479
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCGACCCCTATCCCAGAAGCGGTGGCCGAGGTGAAGTTAAGTGCAGGGGAGGCCGTGGAGCCGGCG
 TCCCTGGGATTTCTGCCCGGGCTGATGGTCAAGGCCCCCTCCCGACCCCTGCCTGCAGGCGCAGCGCAA
 AGAGAAGTCCCGCAACGCGGCGCGCTGGCGGCGGAAAGGAAAACTCGAGTTCTTCGAGCTGGCCAAG
 TTGCTCCCTCTGCCCCGTGCCATCTCCAGTCAGCTGGACAAGGCGTCCATAGTGCCTTACGCTCACAT
 ACCTCCGCTGCGTCGGTTTGGGCTCTGGGAGCGCTCCCTGGGGTTGCGGGCGTCGGGCTCCCGC
 TGGCTCGCCCCAGGTCGCCGAGGCGCTGTGGCTCTGGTGTCTGAAGTCTTCGAGCAACCTAGGTGGA
 CACATCTACAGTCCTTGGATGGCTTCGTGTCGCTTGAACCAGGAAGGAAATTCCTCTACATCTCAG
 AGACAGTGTCCATCTACTTGGTCTCTCACAGGTGGAGCTGACGGGCAGCAGCTCTTCGACTACATCCA
 TCCTGGGACCACTCGGAGGTCTAGAGCAACTTGGACTGAGGGCCGCAACCACTGGCCCCCTACGCCG
 CCTTCTGTGTCTCCACTTCCTCTCGTCGTCTTCTCTCGCTGGTGGACACCCCGAGATTGAAGCCA
 GCCCCACGGAGGCATCCCTGCCTTCCGGTCCAGGAGCGGTCTTCTTTGTCCGAATGAAGTCTACCTT
 CACCAAGAGGGGTCTGAATGTCAAAGCCTCGGGGTACAAGGTCAATCATGTGACAGGACGCTGAGGGCC
 CGAGCGCTGGGTCTTGTAGCCCTTGACACACGCTGCCCCAGCCCCGCTGGCTGAGCTGCCTCTGCATG
 GACACATGATTGTGTTTCGCCTCAGCCTGGGCCTCACCATCCTTGCTTGTGAGAGCAGAGTGAGCGACCA
 TATGGACATGGGGCCCTCAGAACTGGTGGGACGACGTGCTATCAGTTTGTCCACGGACAGGATGTGACC
 AGGATCCGCCAAAGCCATCTGGACCTGCTGGACAAAGGGCAGGTGGTACTGTTACTACCGTTGGCTGC
 AGCGTCGGGGGGCTTCGTGTGGCTCAGTCTGTAGCCACTGTGGCCGGGAACGGGAAGAGTGTGGGA
 GCATCACGTGCTGTGGTCAGTCACGTGCTCAGCCATGCTGAAGGTGGTCAAACACCCCTGGATGCCTTC
 CAGCTTCCAGCCATTGTGTCTCAGGAGGAGTCATCCAGCCAGGCCAGAGCCCAAAAGGAAGAGCCTC
 CAGTTGAGCGAAACAGACTGTGCCTGTGGAACAGGACAAGGACAAGGACCCCTCAGGCCCGAGGCAAACG
 CATCAAAGTGAGGCCAGCCGAAGGAAGTTAGAGGCTCAGAAGACAGTGAGAGGAAGAGGTCTCGGAT
 CCACCGGCTCCGCTCGGCCGAATTCACTTCTGTCTTATAGAGCAGGGGCCCTGAAGCATGATCCAGTAC
 TGCCGTGGGGTTTGAACTCTGGAGACCCCTCGCTGCCCTCCTTCATGCAGGCTTCTGCCACCCGT
 TGTGCGAGGCTGTGTACCCAGGCACCATCCGTTATGGCCCTGCAGAGTTAAGCCTGATGTATCCACAT
 CTGCACAGGCTGGGTCCAGGCCCTCGCTTCCAGAGGCCCTTTACCTACGCTGGGCTGCCCTATCCGG
 GCCCCACAGGTAAGGTGCAGCGAAAGGAGAC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RR210650 representing NM_001107479
 Red=Cloning site Green=Tags(s)


MATPYPRSGGRGEVKCGGGRGAGVPWDFLPGLMVKAPPGPCLQAQRKEKSRNAARWRRGKENLEFFELAK
 LLPLPGAISSQLDKASIVRLSVTYLRLRRFAALGAPPWGLRAVGPPAGLAPGRRGAVALVSEVFEQHLGG
 HILQSLDGFVFALNQEGKFLYISETVSIYLGLSQVELTGSSVFDYIHPGDHSEVLEQLGLRAATTGPPTP
 PSVSSTSSSSSSSLVDPEIEASPTASPAFRVQERSFFVRMKSTLTKRGLNVKASGYKVIHVTGRLRA
 RALGLVALGHTLPPAPLAELPLHGHMIVFRLSLGLTILACESRVDHMDMGPSSELVGRSCYQFVHGQDVT
 RIRQSHLDLLDKGQVVTGYRWLQRAGGFVWLQSVATVAGNGKSAGEHVLWVSHVLSHAEGGQTPDAF
 QLPAAIVSQEESSPGPEPTKEEPPVERKQTPVPEQDKDKDPQARGKRIKVEASPKKEVRGSEDSGEEVSD
 PPAPRPEFTSVIRAGALKHDPVLPWGLTTPGDPSPALLHAGFLPPVVRGLCTPGTIRYGAELSLMYPH
 LHRLGPGPSLPEAFYPTLGLPYPGPTGTRVQRKGD

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001107479

ORF Size: 1785 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001107479.1](#), [NP_001100949.1](#)

RefSeq Size: 1965 bp


RefSeq ORF: 1788 bp

Locus ID: 308387

Cytogenetics: 1q21

MW: 64 kDa

Product images:



Circular map for RR210650