

Product datasheet for **RR210701L3V**

Psmb9 (NM_012708) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Psmb9 (NM_012708) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Psmb9
Synonyms:	Lmp2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012708
ORF Size:	657 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210701).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012708.2 , NP_036840.2
RefSeq Size:	913 bp
RefSeq ORF:	660 bp
Locus ID:	24967
UniProt ID:	P28077
Cytogenetics:	20p12
Gene Summary:	subunit of the proteasome complex, which may play a role in protein catabolism [RGD, Feb 2006]



[View online »](#)