

Product datasheet for **RR210718**

Gpr143 (NM_001106958) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gpr143 (NM_001106958) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Gpr143
Synonyms:	RGD1565799
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR210718 representing NM_001106958 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCCTCCCCGCGCCTGGGAATCTTCTGCTGCCCTTCGTGGGATGCAGCCACACAGCTGGTGCTGACCT
TCCAACCGCGGGTGTCCATGCGCTGTGTCTGGGCAGCGGCCCTCCGCTGGTGCTTGGCCTCCTTCA
GCTCCTAACAGGGCGCCGATCTGTTGGTCACAGGGCGCCTGCGACAACCCAGCAGCCTCAGTCCACATC
CTTCGTGCTGCCACCGCCTGTGATTTGCTTGGCTGCCTGGGAATCGTTATCAGGTCCACAGTGTGGATAG
CCTACCCAGAATTCATTGAAAACATTTCCAATATGAATGGAACAGACATTTGGCCTACTGCTTTCTGTGT
CGGGAGTGCAATGTGGATCCAGCTGTTGTACAGTGCCTGCTTCTGGTGGCTCTTCTGCTATGCAGTTGAT
GTATACTTGGTGATCAGGAGATCAGCAGGACGGAGCACCATCCTGCTGTACCACATCATGGCCTGGGGCC
TGCCTGTGCTGCTCTGTGTGGAAGGTGCAGTCATGCTTTATTACCCTTCTGTGTCCAGGTGTGAGAGAGG
CCTGGACCATGCCATCCCCATTATGTCACCACATACTTGCCACTTATGCTTGTCTTGTGGCAACCCG
ATCCTGTTTACAAGACAGTGATTTTCAGTGGCCTCTTTACTGAAAGGACGAAAAGGTGTTATACAGAGA
ATGAGAGATTGATGGGGCCGTGATCAAGACCCGGTTTTTCAAATAATGCTGGTGTAAATTGCATGTTG
GTTGTCCAATATCATCAATGAATGTCTTTTGTCTACCTTGAATGCAACCAGATACCCATGGAGGCTCT
CTGAAACGCATCCAGAATGCAGCCAGGACCACATGGTTTATTATGGGAATATTGAATCCATCTCAAGGAC
TTCTCTTGTCTCTGGCCTTCTATGGCTGGACAGGATGCAGCCTGGATGTCCATGCTCCAAAGATGGTGAT
TCAGTGGGAAACAATGACTGCCTCGGCTGCTGAGGGCACATATCAGACCCCTGAGGGTCTCTGTGTGCC
CATCAAACCCAGGAAGGTGGTGTGTGTTGGGGGCACACTTCTGATGAGGTGCTGAGTATTTGTCTG
AAGGTTCCGATGCTAGCACTGTTGAAATCCATACTGCAACTGGTCCCACAACATAAAGGAAGTTGACTC
CATTTCCAAGCCAGGGGGATCTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR210718 representing NM_001106958
 Red=Cloning site Green=Tags(s)


MASPRLGIFCCPSWDAATQLVLTFQPRVFHALCLGSGALRLVLGLLQLLTGRRSVGHRAPTTPAASVHI
 LRAATACDLLGCLGIVIRSTVWIAYPEFIENISNMNGTDIWPATFCVGSAMWIQLLYSACFWLFCYAVD
 VYLVIIRRSAGRSTILLYHIMAWGLPVLLCVEGAVMLYPPSVSR CERGLDHAIPHYVTTYLPLMLVLVANP
 ILFHKTVISVASLLKGRKGVYTENERLMGAVIKTRFFKIMLVL IACWLSNIINECLLFYLEMQPDTHGGS
 LKRIQNAARTTWFIMGILNPSQGLLLSLAFYGTGCSLDVHAPKMWIQWETMTASAAEGTYQTPEGSCVP
 HQNPRKVVCVGGHTSDEVL SILSEGSDASTVEIHTATGSHNIKEVDSISQAQGD L

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001106958

ORF Size:	1215 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001106958.1 , NP_001100428.1
RefSeq Size:	1612 bp
RefSeq ORF:	1218 bp
Locus ID:	302619
Cytogenetics:	Xq13
MW:	44.5 kDa