

## Product datasheet for **RR210734L4V**

### Cntn2 (NM\_012884) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cntn2 (NM_012884) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cntn2
Synonyms:	TAG-1; TAG-564; TAG1; Tax
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012884
ORF Size:	3120 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210734).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_012884.1</a> , <a href="#">NP_037016.1</a>
RefSeq Size:	5040 bp
RefSeq ORF:	3123 bp
Locus ID:	25356
UniProt ID:	<a href="#">P22063</a>
Cytogenetics:	13q13
Gene Summary:	a neural cell adhesion molecule that promotes neurite outgrowth [RGD, Feb 2006]



[View online »](#)