

Product datasheet for **RR210736**

Chfr (NM_001009258) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Chfr (NM_001009258) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Chfr
Synonyms:	MGC95316
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RR210736 representing NM_001009258
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAGCCACACGGGAAGACGAGCAGCCGCCACCGGGAGCCCTGGGGAGGCTTCTTCGCCTAGGCG
 CAGAGGAGGACGAACCCAGATCCTTTTGGAAACGCGAGTGGACCATCGGAAGGAGGAGAGGCTGTGA
 CCTCTCTTTTCCAGCAATAAACTGGTCTCTGGAGATCACTGTAACCTTAGAGTGGATGAAAAATCTGGT
 GAGGTGACACTGGAAGACACCAGCACCAATGGAACAGTAATCAACAACTGCGAGTTGTTAAGAAGCAGA
 CTTACCCTTTACAGAGTGGGGATATCATCTATTTGGTGTACAGGAAGAACGAGCCAGAGCACAATGTGGC
 ATACCTCTATGAATCTTTAAATGGCAAACAGAGCTTAACTCAAGAATCCCTTGAAGCCAATAAGAAAAAT
 ATGTTCCATGTGACCAAGAATCCTCAGGTCCAGGGCAGGGTGTGACCCCGAGTTCCACTCTTGTGCAC
 CCATGGCTCAGACATGCTTAGAGGAACACAGCCATCAACATCGACATCAGACCTCCTTCCACGGCCCTC
 TACCTCTTCTACGGAGCCAGAGCTGACCTCTGCAGGGCAAAGCATTCTCTAGCTCTGGACCTGGGACT
 ATAAGCATCTCCCAAAAAGAATGTAGTTCACTTGTGGAAATGGTGAACCTCTAGCCTTTCTTCAGTTT
 TCCAAGACAACGAAGCATCCTTTTCTTTGTTGAAAACCGAAGACACTGAGGAATTGGAGCCTGCCAAAAA
 AAAGATGAAAGGAGATGGGGAGCTTGAACCTCCAGTTATCAGTTGCAGACCAGCGTGAATTGCC
 CAAACCTCAAGTGAAGATGTCAAAGATGCCCTGTGAAACCAGACAAGATGGAGGAGACACTAACCTGTA
 TCATCTGTCAAGACCTTCTGCATGACTGTGTGAGTTTGCAGCCTTGTATGCACACATTTTGTGCGGCTTG
 CTACTCTGGTTGGATGGAGCGTTCATCTCTGTGCCCTACCTGCCGATGTCCAGTGGAGCGGATTTGCAA
 AACCATCCTGAACAACCTTGTGGAAGCATACCTTCTCCAGCACCCAGATAAGAGTCCGAGTGAAGAAG
 ATGTGAGAAGTATGGATGCAAGGAATAAACTCACTCAAGATATGCTGCAACCCAAAGTCAGGAGTCTTT
 CTCTGATGAAGAGGGGAGTTTCCAGAGACCTGTAGAGCTGTCTGATGTTGATAGTGAATCCTCAGACATC
 AGTCAGCCATACATTGTATGCAGACAGTGTCTGAATACAGAAGGCAAGCAGTGCAGTCTTTCCTTGC
 CAGTCCCAGAGAGTGAAGTGGGAGCCACACAGGCCCTTGGTGGGAAGCACCTTCAACATCTGCCAGCTT
 GACAACAGCCCCGGTTATACATGCCCTTTCAGGGAAGCCATGCCATATGCACCTGCTGCTTCCAGCCT
 ATGCCGGACCGGAGAGCTGAGCATGAGCAGGATCCCGCATCGCCCTCAACAGTGTGCAGTGTGCCTGC
 AGCCCTTCTGCCACCTGTACTGGGGCTGCACGCGAACTGGCTGCTTTGGCTGCTGGCCCCATTCTGTGA
 GCTCAACCTGGGGACAAGTGTGGATGGAGTGTGAACAATAACAATTATGAGTCGGATATCCTGAAG
 AATTACCTGACAACCGGGTCTGACATGGAAAAATATGTTGACAGAGAGTCTCCTGGCTCTACAGCGGG
 GTGATTTATACTGTCTGATTACAGAATCACTGGAACACCGTGTGTGTTACTGCTGTGGTCTCGGTAG
 CTTCCGAGAGCTGACCTATCAGTATCGTCAGAATTCCTGCTTCAGAGTTGCCAGTGACCGTAACATCC
 CGTCTGATTGCTACTGGGGCCGCAACTGTCGCACTCAGGTGAAGGCTCACCATGCAATGAAATCAATC
 ACATCTGTGAGCAACAAGGTTCAAGAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR210736 representing NM_001009258
Red=Cloning site Green=Tags(s)

```
MEPHGEDEQPPPREPWGRLRLGAEDEPQILLWKREWTIGRRRGCDLFFPSNKLVS GDHCKLRVDEKSG
EVTLEDSTNGTVINKLRVVKQTYPLQSGDIIYLVYRKNEPEHNVAYLYESLNGKQSLTQESLEANKEN
MFHVTKESSGPGQDDPQVPLLSPMAQTCLEEQPSTSTSDLLPTASTSSTEPELTSAGQKHSSSSGPGT
ISISPKECSSLVANGELSSLSSVFQDNEASFLLKTEDTEELEPAKKMKMGDGLDLNLQLSVADQRGIA
QTSSDVKDASVKPKMEETLTCIICQDLLHDCVSLQPCMHFCAACYSGWMERSSLCPTCRCPVERICK
NHILNMLVEAYLLQHPDKSRSEEDVRSMDARNKITQDMLQPKVRRSFDEEGSDELLELSDVDSESSDI
SQPIYIVCRQCPEYRRQAVQSLPCVPPESELGATQALGGEAPSTSASLTTAPGYTCPLQGSHAICTCCFQP
MPDRRAEHEQDSRIAPQQCAVCLQPFCHLYWGCTRTGCFGCLAPFCELNLGDKCLDGLVNNNNYESDILK
NYLTTTRGLTWKNMLTESLLALQRGVFIILSDYRITGNTVLCCYCGLRSFRELTQYRQNPASELPVTVTS
RPDCYWGRNCRTQVKAHAMKFNHICEQTRFKN
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001009258

ORF Size: 1989 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001009258.1](#), [NP_001009258.1](#)

RefSeq Size: 3185 bp


RefSeq ORF: 1992 bp

Locus ID: 288734

Cytogenetics: 12q16

MW: 74 kDa

Product images:



Circular map for RR210736