

Product datasheet for **RR210746L4V**

Psmc3ip (NM_134458) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Psmc3ip (NM_134458) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Psmc3ip
Synonyms:	Gt198; Tbpip
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_134458
ORF Size:	651 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210746).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_134458.2 , NP_604453.2
RefSeq Size:	936 bp
RefSeq ORF:	654 bp
Locus ID:	140938
UniProt ID:	Q91ZY6
Cytogenetics:	10q31
Gene Summary:	interacts with nuclear receptors, estrogen receptor alpha and beta, thyroid hormone receptor beta1, androgen receptor, glucocorticoid receptor, and progesterone receptor [RGD, Feb 2006]



[View online »](#)