

Product datasheet for **RR210763L3V**

Pdha2 (NM_053994) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pdha2 (NM_053994) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pdha2
Synonyms:	Pdhal
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053994
ORF Size:	1173 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210763).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053994.2 , NP_446446.1
RefSeq Size:	1574 bp
RefSeq ORF:	1176 bp
Locus ID:	117098
UniProt ID:	Q06437
Cytogenetics:	2q44
Gene Summary:	mouse homolog codes for one of two E1 alpha subunits for pyruvate dehydrogenase (Pdh) [RGD, Feb 2006]


[View online »](#)