

Product datasheet for **RR210777L3V**

Zdhhc5 (NM_001039338) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Zdhhc5 (NM_001039338) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Zdhhc5
Synonyms:	DHHC-5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001039338
ORF Size:	2145 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210777).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001039338.1 , NP_001034427.1
RefSeq Size:	2697 bp
RefSeq ORF:	2148 bp
Locus ID:	362156
UniProt ID:	Q2THW7
Cytogenetics:	3q24
Gene Summary:	Palmitoyl acyltransferase for the G-protein coupled receptor SSTR5 and for FLOT2. [UniProtKB/Swiss-Prot Function]


[View online »](#)