

Product datasheet for **RR210796L4V**

Kcnip3 (NM_032462) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Kcnip3 (NM_032462) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Kcnip3
Synonyms:	Csen; Dream; rKChIP3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_032462
ORF Size:	768 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210796).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032462.2 , NP_115851.1
RefSeq Size:	818 bp
RefSeq ORF:	771 bp
Locus ID:	65199
UniProt ID:	Q9JIM47
Cytogenetics:	3q36
Gene Summary:	member of the voltage-gated K ⁺ channel-interacting proteins (KCNIPs) family; plays a role in the intracellular trafficking of K ⁺ channel alpha subunit Kv4.2 (Kcna4); modifies Kv4 potassium channel current kinetics [RGD, Feb 2006]



[View online »](#)