

## Product datasheet for **RR210810**

### **Gjb3 (NM\_019240) Rat Tagged ORF Clone**

#### Product data:

Product Type: Expression Plasmids  
Product Name: Gjb3 (NM\_019240) Rat Tagged ORF Clone  
Tag: Myc-DDK  
Symbol: Gjb3  
Vector: pCMV6-Entry (PS100001)  
E. coli Selection: Kanamycin (25 ug/mL)  
Cell Selection: Neomycin  
ORF Nucleotide Sequence: >RR210810 representing NM\_019240  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGGATTGGAAGAAGCTTCAGGACCTGTTGAGCGGTGTGAACCACTACTCCACCGCATTCCGGGCGCATCT  
GGCTGTCGGTAGTGTTCGTCTCCGGGTGCTGGTGTATGTGGTGCAGAGCGTGTGTGGGGCGACGA  
GCAAAAAGACTTTGACTGTAACACCAGGCAGCCGGTTGTACCAACGTGTGCTATGACAACCTTCTCCCC  
ATCTCCAACATCCGCCTCTGGGCCCTGCAGCTCATCTTCGTACAGTGTCCCTCTATGCTGGTCATCCTGC  
ACGTAGCCTACCGCGAGGAGCGGAAACGGAAGCATCGCCAGAAGCACGGGAGCATTGCGCCAAACTGTA  
CAGCCACCCGGGCAAGAAGCACGGCGCCTGTGGTGGACCTACCTGTTCAAGTCTCATCTTCAAGCTCATC  
ATTGAATTGGTCTTCTGTATGTTCTACACAGCTCTGGCATGGCTTACCATGCCGCGTCTGGTACAGT  
GCGCCAGCGTGGTACCTTGCCCCAACCCGTGGATTGCTACATCGCTCGGCCACGGAGAAGAAAGTCTT  
TACCTACTTCATGGTAGGCGCTGCTGCCGTCTGCATTATCTCACCATCTGTGAGATCTGCTACCTCATC  
TTCCACAGGATTATGCGAGGCCGTGAGCAAGGACAAATCGACGAAGAGCATCAGCTCCCCGAAGTCTCCA  
GCCGGCCTCCACCTGTCGCTGTACCACAAGCTGCTGGAGAGTGGTGTCTGGAAGCAGTACCAGCCGA  
TGACAAGCTGCAGGCTTACGCGCTAGCCTGACCCCAT

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RR210810 representing NM\_019240  
Red=Cloning site Green=Tags(s)

MDWKKLQDLLSGVNQYSTAFGRWL SVVFVFRVLVYVVAERVWGDEQKDFDCNTRQPGCTNVCYDNFFP  
 ISNIRLWALQLIFVTCPSMLVILHVAYREERERKHRQKHGEHCAKLYSHPGKKHGGGLWWTYLFSLIFKLI  
 IELVFLYVLHTLWHGFTMPRLVQCASVVPNPNTVDCYIARPTKVKFTYFMVGASAVCIILTICEICYLI  
 FHRIMRGLSKDKSTKSISSPKSSSRASTCRCHKLLESGDLEAVPADDKLQASAPSLTPI


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_019240

**ORF Size:** 810 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_019240.1</a> , <a href="#">NP_062113.1</a>
<b>RefSeq Size:</b>	813 bp
<b>RefSeq ORF:</b>	813 bp
<b>Locus ID:</b>	29585
<b>UniProt ID:</b>	<a href="#">P25305</a>
<b>Cytogenetics:</b>	5q36
<b>MW:</b>	31 kDa
<b>Gene Summary:</b>	may play a role in hearing [RGD, Feb 2006]