

Product datasheet for **RR210819**

Lrriq3 (NM_001024307) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lrriq3 (NM_001024307) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Lrriq3
Synonyms:	Lrrc44
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RR210819 representing NM_001024307
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTTTCATGGAACAATCACAAGGAGCTTACCAGCCATGAGGAGTGGAGTCATTATAATGAAAATATAG
 TAGAAGATCAGAAAGACTTTGTTTTGTGAAGTATAATGGCCTTCACCTGAAGTCCATGGAGAACCTGCA
 GTCCTGCATCTCTGAGGGTCTGCATCTTCTCTAACAACCTTTGTCACTGACATTCAGCCACTGCAGGGC
 TGCAAAAAGCTCATCAAACCTTGATCTCCATGGAACCAGATAAAGACTCTTCAGATAGGACATTTTGGA
 ATGGACTCAAGAACCCTCAAACCTCTATCTGCATGACAATGGGTTTGCCAACTGAAGAACATCTGTGT
 ATTGTCTGGTTGTGTAACCTCGTAGGCCTGACCATGTTTGACTGCCAGTGAGCCTGAAGAAAGGCTAC
 AGGCATGTTCTTGTAAGTATCTGGCCTCTAAAAGCTCTGGACCACCACGTGATATCCGATGAAGAAA
 TCATTCAGAACTGGCATCTTCAGAGAGATTTAAACGTTAGCCAAAGTCTATTCTTTAACATTTGCC
 AGCACCAGATAAAGGGCACAACCTATGAAGATGAAATAAAAAATATCAAGCATATAATTTCAAGGATCAAT
 GAAATCCTGGCTCATAACTCACCGGTATTGATCATCCAGAGGTGGATCCGTGGTTTCATAGTTAGGAAAC
 ACTTAAGCCCCTATTTTACACGTAAGGATCACGACAGAATGATTAGGTAAGTGAACAAAAATTTGAT
 CTGCATAGGTAGAAATAACGAAGACAAATCTTGAGGATATCTTCTTATCAAGCCTGAATCTAATATA
 AAAGGGAAGGTCGACACTGGAACAAATGAGATATTCTCTGATGACTTTAAATATTCTACTGAGTACA
 GAAAACACATTTCTGTCTTTCATATGAGTTGAAAACAAAAGATATTGCTGGAAATTCAAAACAACCAAG
 ACATCACATTCAGAAAGGTCAGAACGAGATGAAATCTGATTCTGAGGATGAGGAAGTGGACCAAGTTT
 AGGATATCTGCGATGAAGATCCCTCTATGTCCTTCCAGATCACTCAAGTACCGGGCAATGCTGAAGGAGA
 TAAAATGGGACTACTTCCCGGAGTATCTCCAGCCACTCCCTGCCACTCGTCAAAAGCCACAAGTTAAAAG
 AGAGACTCATGAGGAGCTGAAGGAAAGGAGAAATTTTCTAGCCTCACAGAGGGCTGGCATGGATCTCCAC
 TTGTTCAATGATATTGATAAATATTACTCAGAGCAAAGGCAGCAAGAGGAAGAAGCAAGAAAATTTGCAG
 CCACAGCTGCAGCACAAGTTGCCAGGAAAGAGCTAGGCTCAACACTAGTGAAAATTTAAGGAAGAAGAT
 TTCTATGGCAAGAAAATAATGCAAAAGGATAATGAGACAATCCAGAAAGGCTTACGGCACATCTGGAAG
 GAGAGACTCACTACCTGGAAAAGGTTAGAGAGCGGAAATCCATGTTTCTAGCTGAAAAAACTGAATG
 CTGCAGACCAGTCTTTAGTTCTTTCCCTCAACAATGAGAGGTCCATTCTGCTCAGAGGAATCATTCAAGT
 TGACAAACTGAAGCGAGACTTGCTGAAATGAAAGCAAAGCACTTCGATGTTATTGAGAAACGGGAAAAA
 CGAAAATACAAGCAGAATTTGTTGACGGAAATGAAGAAGCTAAGGCGGAGGAAATCCATAAAGGCACT
 GTGAAGAGAGATTTGTCAATTGATACGCTCATCTTTCAAAAAGGCTGTGAAGACTTGAGGAGGCGAAAGC
 AAAAGTTGAATTTATTAACACATATTACAGTCCAAATCTCACAAAAATATCATAATAAAACAATGCTC
 AATCACATC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR210819 representing NM_001024307
 Red=Cloning site Green=Tags(s)


MFHGTITKELTSHEEWSHYNENIVEDQKDFVFKYNGHLKSMENLQSCISLRVCIFSNNFVTDIQPLQG
 CKKLIKLDLHGNQIKTLPDRTFWNGLKNLKL LYLHDNGFAKLKNICVLSGCVNLVGLTMFDCPVSLKKG
 RHVLVNSIWPLKALDHHVISDEEIIQNWHLPERFKTFSQSLFFNICPAPIKGTTYEDEIKNIKHII SRIN
 EILAHNSPVLIIQRWIRGFIVRKHLSPYFTRKRHHDRMIRVLETKLICIGRNNEKFLIEDIFLIKPE
 SNI KGKVAHWKQMRYSRDDFKYSTEYRKHISCLSYELKTKDIAGNSKQPRHHIQKGQNMKSDSEDEEVD
 TKF RISAMKIPLCPSRSLKYRAMLKEIKWDYFPEYLQPLPATRQKPQVKRETHEELKERNFLASQ
 RAGMDLH LFNIDIDKYSEQRQQEEEEKFAATAAAQVAQERARLNTSENLRKKISMARKLMQK
 DNETIQKGLRHIWK ERLNYLEKVRERKSMFLAEKKLNAADQSLVLSLNNERSILLRGIQV
 DKLKRDLSMKAKHFDVIEKREK RKYKQNLLETMKLRAEEIHKRHCEERFVIDTLIFQKGCER
 LEEAKAKVEFINTYYTSKSHKKYHNKTM
 L NHI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001024307

ORF Size: 1899 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001024307.1](#), [NP_001019478.1](#)

RefSeq Size: 2417 bp

RefSeq ORF: 1902 bp


Locus ID: 499732

UniProt ID: [Q6AYL8](#)

Cytogenetics: 2q45

MW: 75.1 kDa

Product images:



Circular map for RR210819