

Product datasheet for **RR210862**

Rab8a (NM_053998) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rab8a (NM_053998) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Rab8a
Synonyms: Mel; Rab8
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR210862 representing NM_053998
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCGAAGACCTACGATTACCTGTTCAAGCTGCTGCTGATCGGGGACTCGGGGGTAGGGAAGACCTGTG
TCCTGTTCCGCTTCTCCGAGGACGCCTTCAATTCACGTTTCATCTCTACCATAGGAATTGACTTTAAAT
TAGGACCATAGAGCTCGATGGCAAGAGGATTAAGTGCAGATATGGGACACAGCCGCGCCAGGAGCGGTTT
CGGACAATCACGACAGCCTACTACAGGGGTGCCATGGGTATCATGCTGGTCTACGACATTACCAATGAGA
AGTCCCTTCGACAACATCCGGAATTGGATTCCGGAACATTGAGGAGCATGCCTCTGCAGATGTGGAGAAGAT
GATACTGGGGAATAAGTGTGATGTGAACGACAAGAGGCAGGTGTCCAAGGAACGGGGAGAAAAGCTGGCA
CTCGACTATGGGATCAAGTTCATGGAGACCAGCGCAAGGCCAACATCAATGTGGAGAACGCATTTTCA
CTCTCGCCAGGGATATCAAAGCAAAAATGGACAAAAATTTGGAAGGGAACAGCCACAGGGGAGCAGCCA
TGGAGTCAAGATCACAGTGGAGCAGCAGAAGAGGACCAGCTTCTCCGGTGCAGTCTCTCTG

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR210862 representing NM_053998
Red=Cloning site Green=Tags(s)


MAKTYDYLFKLLLIGDSVGVKTCVLFVRFSEDAFNSTFISTIGIDFKIRTIELDGKRIKLQIWDTAGQERF
RTITTAYYRGAMGIMLVYDITNEKSFNIRNWIRNIEEHASADVEKMILGNKCDVNDKRQVSKERGEKLA
LDYGIKFMETSAKANINVENAFFTLARDIKAKMDKKLEGNSPQGS SHGVKITVEQQRTSFFRCSLL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



Cloning Scheme:

Plasmid Map:


ACCN: NM_053998

ORF Size: 621 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_053998.2](#), [NP_446450.2](#)

RefSeq Size: 1768 bp

RefSeq ORF: 624 bp

Locus ID: 117103

UniProt ID: [P35280](#)

Cytogenetics: 16p14

MW: 23.7 kDa

Gene Summary: member of low molecular weight (LMW) GTP-binding protein, rab family; may function in regulating neurotransmitter release [RGD, Feb 2006]