

Product datasheet for **RR210864**

Slc27a1 (NM_053580) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Slc27a1 (NM_053580) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Slc27a1
Synonyms:	FATP
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RR210864 representing NM_053580
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGCGGACTCCGGGAGCAGGAACAGCCTCTGTGGCCTCATTGGGGCTGCTTTGGCTTCTGGGACTTCCGT
 GGACCTGGAGCGCGCGCGCGCTTCGGTGTGTACGTGGGTAGCGGTGGCTGGCGCTTCTGCGTATCGT
 CTGCAAGACGGCGAGGCGAGACCTCTTTGGCCTCTCTGTTCTGATCCGCGTGCGGCTAGAGCTACGACGA
 CACCGGCGAGCAGGAGACACGATCCACGCATCTTCCAGGCCGTGGCCAGCGACAGCCGGAGCGCTGG
 CGCTGGTAGATGCGAGTAGCGGTATCTGCTGGACCTTCGCACAGCTAGACACCTACTCCAATGCTGTGGC
 CAATCTGTTCTCCAGCTGGGCTTTGCGCCAGGCGATGTGGTGGCTGTGTTCTGGAAGGCCGGCCGAG
 TTCGTTGGGACTGTGGCTGGCCAGGCCGGTGTAGTGGCTGCGCTTCTCAATGTCAACCTGAGGC
 GGGAGCCCCTTGCTTCTGCTTGGGCACATCAGCTGCCAAGGCCCTCATTTATGGCGGGGAGATGGCAGC
 GCGGTTGGCGGAGGTGAGTGAGCAGCTGGGGAAGAGCCTGCTCAAGTTCTGCTCTGGAGATCTGGGGCT
 GAGAGCGTCTGCTGACACGCAGCTTCTGGACCCATGCTTGTGAGGGGCCACCACACCCCTGGCAC
 AGGCCCAAGGCAAGGCGATGGATGATCGCTATTTTACATCTATACTTCTGGGACCACCGGACTTCTCAA
 GGCGGCCATTGTGGTGCACAGCAGGTAACCGCATTGCAGCCTTCGGCCACCATTCTACAGCATGCGG
 GCCAACGATGTGCTCTATGACTGCCTACCTCTCTACCACTCAGCAGGGAACATCATGGGCGTGGGACAGT
 GTATCATCTACGGGTTAACGGTGGTACTGCGCAAGAAGTTCTCCGCCAGCCGCTTCTGGGACGACTGTGT
 CAAATAAATGACCGGTAGTGCAGTACATCGGTAAAATATGCCGCTACCTGCTAAGGCAGCCGGTTCCG
 GATGTAGAGCGCGGCACCACGTGCGCCTGGCCGTGGGTAACGGACTGCGGCCAGCCATCTGGGAGGAGT
 TCACGCAGCGTTTCGGTGTGCGACAGATTGGCGAGTTCTACGGCGCCACCAGTGAACCTGCAGCATTGC
 CAACATGGACGGCAAGGTGCGCTCCTGCGGCTTCAACAGCCGTATCCTCACGCATGTGTACCCCATCCGT
 CTGGTCAAGGTCAACGAGGACACGATGGAGCCACTGAGGGACTCCCAAGGCCTCTGCATCCCGTGGCAGC
 CCGGGGAACCTGGCCTTCTCGTGGGCCAGATCAACCAGCAAGACCCTCTGCGGCGCTTCGATGGCTATGT
 TAGTGACAGCGCCACCAACAAGAAGATTGCCACAGCGTGTTCGAAAAGGGGACAGCGCTACCTTTCA
 GGTGACGTGCTAGTGTGGACGAGCTGGGTACATGTACTTCCGTGACCGCAGCGGGGATACCTTCCGAT
 GGCGCGGCGAGAACGATCCACCACGGAGGTGGAAGCCGTGCTGAGCCGCTGTTGGGCCAGACGGACGT
 GGCTGTGTATGGAGTGGCTGTGCCAGGAGTGGAGGGAAAGCCGGCATGGCGGCCATTGCAGACCCAC
 AACCAGCTGGACCCTAATCAATGTACCAGGAATTGCAGAAGGTTCTTGCATCCTATGCCAGCCCATCT
 TCCTGCGTCTTCTGCCCAAGTGATACAACAGGCACCTTCAAGATCCAGAAGACCCGGCTACAGCGTGA
 AGGCTTTGACCCCGCCAGACCTCAGACCGGCTCTTCTTCTAGACCTGAAACAGGGACGCTACCTACCC
 CTGGATGAGAGAGTCCATGCCCGCATCTGCGCAGGCGACTTCTCACTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR210864 representing NM_053580
 Red=Cloning site Green=Tags(s)


MRTPGAGTASVASLGLLWLLGLPWTWSAAAAFGVYVGGSGWRFLRIVCKTARRDLFGLSVLIRVRELR
 HRRAGDTIPRIFQAVAQRQPERLALVDASSGICWTF AQLDTYSNAVANLFLQLGFAPGDVVAVFLEGRPE
 FVGLWLGLAKAGVVAALLNVNLRREPLAFCLGTSAAKAL IYGGEMAAVAEVSEQLGKSLKFCSDGLGP
 ESVL PDTQLLDPMLAEAPTTPLAQAPGKGMDDRLFYIYTS GTTGLPKAAI VVHSRYRIA AAFGHHSYSMR
 ANDVLYDCLPLYHSAGNIMGVQCIIYGLTVVLRKKFSASRFWDDCVKYNCTVVQYIGEICRYLLRQPVR
 DVERRHHVRLAVGNLRPAIWEFFTQRFVGRQIGEFYGATECNC SIANMDGKVGSGCFNSRIL THVYPIR
 LVKVNEDTMEPLRDSQGLCIPCQPGEPGLLVGQINQDPLRRFDGYVSDSATNKKIAHSVFRKGD SAYLS
 GDVLMDELGYMYFRDRSGDTFRWRGENVSTTEVEAVLSRL LGQTDVAVYGVAVPGVEGKAGMAIADPH
 NQLDPNSMYQELQKVLASYAQPIFLRLLPQVDTTGTFKIQKTRLQREGFDPRTSDRLF FLDLKQGRYLP
 LDERVHARICAGDFSL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_053580

ORF Size: 1938 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_053580.2](#), [NP_446032.2](#)

RefSeq Size: 3210 bp
RefSeq ORF: 1941 bp
Locus ID: 94172
Cytogenetics: 16p14
MW: 71.3 kDa
Gene Summary: involved in myocardial fatty acid uptake [RGD, Feb 2006]