

Product datasheet for **RR210872L4V**

Pgr (NM_022847) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pgr (NM_022847) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pgr
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_022847
ORF Size:	2769 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210872).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022847.1 , NP_074038.1
RefSeq Size:	3056 bp
RefSeq ORF:	2772 bp
Locus ID:	25154
UniProt ID:	Q63449
Cytogenetics:	8q11
Gene Summary:	receptor for progesterone; regulates glutamic acid decarboxylase expression in the hypothalamus during proestrus; may mediate production of the GnRH and LH surge [RGD, Feb 2006]



[View online »](#)