

Product datasheet for **RR210919L4V**

Paics (NM_080910) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Paics (NM_080910) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Paics
Synonyms:	Ade2h1; Airc
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_080910
ORF Size:	1275 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210919).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_080910.3 , NP_543186.1
RefSeq Size:	1994 bp
RefSeq ORF:	1278 bp
Locus ID:	140946
UniProt ID:	P51583
Cytogenetics:	14p11
Gene Summary:	bifunctional enzyme that catalyzes the 6th and 7th step in de novo purine nucleotide synthesis [RGD, Feb 2006]



[View online »](#)