

Product datasheet for **RR210920L4V**

Dusp1 (NM_053769) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dusp1 (NM_053769) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Dusp1
Synonyms:	3CH134; CL100; MKP-1; Mkp1; Ptpn16
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053769
ORF Size:	1101 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210920).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053769.3 , NP_446221.2
RefSeq Size:	1940 bp
RefSeq ORF:	1104 bp
Locus ID:	114856
UniProt ID:	Q64623
Cytogenetics:	10q12
Gene Summary:	This gene encodes a protein that catalyzes the dephosphorylation and inactivation of MAP kinase, and may be involved in insulin mediated signalling. [provided by RefSeq, Jul 2008]



[View online »](#)