

## Product datasheet for **RR210948L3V**

### Zfp111 (NM\_133323) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Zfp111 (NM_133323) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Zfp111
Synonyms:	rKr2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_133323
ORF Size:	2232 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210948).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_133323.1</a> , <a href="#">NP_579857.1</a>
RefSeq Size:	2817 bp
RefSeq ORF:	2235 bp
Locus ID:	170849
Cytogenetics:	1q21
Gene Summary:	has transcriptional repressor activity; may act as a transcriptional repressor in oligodendrocytes of the CNS; may play a role in myelination [RGD, Feb 2006]



[View online »](#)