

Product datasheet for **RR210962**

Tmbim6 (NM_019381) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tmbim6 (NM_019381) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Tmbim6
Synonyms:	Ab1-011; Ac1-149; Cc1-27; Tegt
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR210962 representing NM_019381 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAATATATTTGATCGGAAGATCAACTTTGATGCCCTCTAAAATTTCCACATAACTCCCTCCACAC
AGCAGCACCTAAAGAAGGTCTATGCCAGTTTTGCACTGTGCATGTTTGTGGCAGCAGAGGGCCTATGT
CCATGTGGTCACACGTTTCATCCAGGCTGGCCTGCTCTGCCCTGGGCGCCCTGGCCTTGATGATTTC
CTGATGGCCACACCTCACAGCCATGAGACGGAGCAGAAGAGGCTGGGACTGCTCGCTGGCTTCGCCTCC
TTACAGGAGTTGGCCTGGGACCTGCCCTGGAGCTGTGCATTGCCATCAACCCAGCATCCTCCCCACGGC
CTTCATGGGCACGGCCATGATCTTCACCTGCTTCAGCCTGAGTGCCTCTACGCCAGGCGCCGGAGTTAC
CTCTTTTTGGGAGGTATCTTGATGTGAGCCATGAGCCTCATGTTTCGTGTCCTCTCTGGGAACTTTTCT
TTGGATCCATTTGGCTGTTCCAGGCAAACCTGTACATGGGGCTGCTGGTCATGTGCGGCTTTGTCCTCT
CGACACTCAGCTCATTATTGAAAAGGCTGAACACGGAGACAAGGATTACATCTGGCACTGCATTGACCTC
TTCTTGGACTTCGTTACACTCTTCAGGAAGCTCATGCTGATCTTGGCCTCAATGAGAAGGACAAAAAGA
AAGAGAAGAAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR210962 representing NM_019381
 Red=Cloning site Green=Tags(s)

MNIFDRKINFDALLKFSHITPSTQQHLKKVYASFALCMFVAAAGAYVHVVTRFIQAGLLSALGALALMIC
 LMATPHSHETEQRLLAGFAFLTGVGLGPALELCIAINPSILPTAFMGTAMIFTCFSLSALYARRRSY
 LFLGGILMSAMSLMFVSSLGNLFFGSIWLFQANLYMGLLVCMGFVLFDTQLIEKAEHGDKDYIWHCIDL
 FLDFVTLFRKLMILAFNEKDKKKEKK


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_019381
ORF Size: 711 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_019381.2 , NP_062254.2
RefSeq Size:	898 bp
RefSeq ORF:	714 bp
Locus ID:	24822
UniProt ID:	P55062
Cytogenetics:	7q36
MW:	26.5 kDa
Gene Summary:	Suppressor of apoptosis [RGD, Feb 2006]