

Product datasheet for **RR210969**

Vsig2 (NM_001106812) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Vsig2 (NM_001106812) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Vsig2
Synonyms: RGD1561274
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR210969 representing NM_001106812
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCCTGGCCCTGGTGGGGCTTTCCTTTGTGGTCACTTGCTGGGCTTCTGTGCTTGAGCGGGCTAG
CTGTGGAGGTGACTGTACCCGCTGAGCCAGTGAGCATGCCAATGGCAAGACGGCCGAGCTGAGCTGCAG
TTACAGCACGTCGGTGGGCGATAACTTTGCCTTGGAGTGGAGCTTGTGCAACCTGGGAAGCCATCTCT
GCATCTGCGCCTATTCTGTACTTCACCAACGGCCACCTGTATCCCACTGGCTCTAAAGCAGATCGGGCCA
GCCTTCTCCACAACCCGCCACAGGAGGACTGGCTACCCTGAAACTGACTGACCTCCGCCCTCAGACAC
TGGAACCTACCTCTGCAATGTTAACAACCCACCAGACTTCTACACCAACGGTTTGGGGTTAATCAACCTT
ACCGTGTCTGGTGCCCCAGTCACCCCTTATGCAGCCAGAGTGGGCAAACCTCTGTGGGAGGTTCTGCTG
CACTGAGATGCAGCTCTTCTGAGGGAGCACCGAAGCCAGTGTACAACGGTACGTCTTGGATCTTCTCC
TACACCTCCTCCGGTAGCATGGTTCAAGATGAAGTGTCTGGCCGGCTCATTCTACCAATCTCTCCCTG
ACCTCCTCTGGTACCTACCGATGTGTAGCCAGCAACAGATGGCAGTGCCTCCTGTGAGCTGAACCTCT
CTGTGACTGACTCATCTGAAGCCGAGTGGCTGGGACTCTGATTGGGGTGTCTCTGGGTGTCTGTTGCT
GTCGGTTGCTGCATTCTGCCTGATCAGGTTCCAGAAAGACAGAAAAAGGAGCCCAAGGAGACATATGGG
GGTAGTGACCTTCGGGAGGATGCCACCGCTCCTGGGGTTTCTGAGCAAGCTTCTGTGAGAGCAGACTGTA
GCAAAGGGCTCCTGGAGAAGTCTCCAGCTGCCAGCACTGTGACCACCACCAATCCAAACTCTATGGT
TGTC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR210969 representing NM_001106812
 Red=Cloning site Green=Tags(s)

MAWPLVGAFLCGHLLGFLCLSGLAVEVTVPAEPVSMMPMGKTAELSCSYSTSVGDNFALEWSFVQPGKPIS
 ASAPILYFTNGHLYPTGSKADRASLLHNPPTGGLATLKL TDLRPSDTGTYL CNVNNPPDFYTNGLGLINL
 TVLVPPSHPLCSQSGQTSVGGSAALRCCSSSEGAPKPVYVWVRLGSSPTPPPGSMVQDEVSGRLILTNLSL
 TSSGTYRCVASNQMGASCELNLSVTDSSSEGRVAGTLIGVLLGVLLLLSVAAFCLIRFQKDRKKEPKETYG
 GSDLRE DATAPGVSEQASVRADCSKGLLEKSPAASTVTTT KSKLSMVV


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001106812

ORF Size: 984 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001106812.1 , NP_001100282.1
RefSeq Size:	1175 bp
RefSeq ORF:	987 bp
Locus ID:	300520
Cytogenetics:	8q21
MW:	34.2 kDa