

Product datasheet for **RR210978L4V**

Khdrbs1 (NM_130405) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Khdrbs1 (NM_130405) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Khdrbs1
Synonyms:	P62; Sam68
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_130405
ORF Size:	1329 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR210978).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_130405.1 , NP_569089.1
RefSeq Size:	2655 bp
RefSeq ORF:	1332 bp
Locus ID:	117268
UniProt ID:	Q91V33
Cytogenetics:	5q36
Gene Summary:	RNA binding protein and acts as a substrate of the Insulin receptor and may have a role as a docking protein in Insulin receptor signaling [RGD, Feb 2006]



[View online »](#)