

Product datasheet for **RR210981**

Egfl7 (NM_139104) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Egfl7 (NM_139104) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Egfl7
Synonyms: MGC93163
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR210981 representing NM_139104
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCCGAGGGTTCTGGAACTCCCTGAATGGGTGACAGACAGGCCTCTCTACCCACAGCAATATGCC
AGCCACCGTGTGGGAATGAAGGGAGTTGCATCCGCCAGGCCGCTGCCGCTGCCCTGTGGATGGCAGGG
AGATACTTGCCAGATAGATGTTGATGAATGCAGTACAGGAGAGGCCAGGTGTCCCAGCGCTGTGTCAAT
ACTGTGGGAAGTTACTGGTGCCAGTGTGGGAGGGGCAGAGCCATCTGCAGATGGGGTGTCTGCCTGC
CTAAGGAGGGGCCCTCCCGGTGGCACCAGTCCACGCCAGGAGTGGACAGCGTGGTGAGAGAGRAGGT
GTACAAGCTGCAGGCCGGGTTGATGTATTAGAACAGAACTGCAGCTGGTGTGGCCCCACTGCACAGC
CTGGCCTCTCGGTCCCGGAGCATGGGCTGCAAGATCCTGGCAGCCTGCTGGCCCCACTCCTCCAGCAGC
TGGACCGATTGACTCACTGAGCGAGCAGGTCTCCTTCTGGAGGAACAGCTGGGGTCTGCTCTTGCAA
AAAAGACCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR210981 representing NM_139104
Red=Cloning site Green=Tags(s)

MPQSGNSPEWVTRPLLPTAICQPPCGNEGSCIRPGRRCRPVWQGDTCQIDVDECSTGEARCPQRCVN
TVGSYWCQWEGQSPSADGVLCLPKEGSPVAPSPTPGVDSVVREXVYKQARVDVLEQKLQLVLAFLHS
LASRSPHEGLQDPGSLLAHSFQQLDRIDSLSEQVSFLEEQLGSCSCKDL


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_139104

ORF Size: 570 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_139104.1](#), [NP_620804.1](#)

RefSeq Size: 1888 bp

RefSeq ORF: 573 bp

Locus ID: 245963

UniProt ID: [Q6AZ60](#)

Cytogenetics: 3p13

MW: 20.5 kDa

Gene Summary: Regulates vascular tubulogenesis in vivo. Inhibits platelet-derived growth factor (PDGF)-BB-induced smooth muscle cell migration and promotes endothelial cell adhesion to the extracellular matrix and angiogenesis (By similarity).[UniProtKB/Swiss-Prot Function]