

Product datasheet for **RR211014L4V**

Ednrb (NM_017333) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ednrb (NM_017333) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ednrb
Synonyms:	Ednra; Etb
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017333
ORF Size:	1326 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211014).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017333.1 , NP_059029.1
RefSeq Size:	2018 bp
RefSeq ORF:	1329 bp
Locus ID:	50672
UniProt ID:	P21451
Cytogenetics:	15q22
Gene Summary:	Non-specific receptor for endothelin 1, 2, and 3. Mediates its action by association with G proteins that activate a phosphatidylinositol-calcium second messenger system. [UniProtKB/Swiss-Prot Function]



[View online »](#)