

## Product datasheet for **RR211021L4V**

### **Rgs4 (NM\_017214) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Rgs4 (NM_017214) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Rgs4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017214
ORF Size:	615 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211021).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_017214.1</a> , <a href="#">NP_058910.1</a>
RefSeq Size:	2919 bp
RefSeq ORF:	618 bp
Locus ID:	29480
UniProt ID:	<a href="#">P49799</a>
Cytogenetics:	13q24
Gene Summary:	regulates G-protein coupled receptor signaling by accelerating the inactivation of active G-proteins; may regulate dopamine D2 receptor signaling [RGD, Feb 2006]



[View online »](#)