

## Product datasheet for **RR211023**

### Il1b (NM\_031512) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Il1b (NM\_031512) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Il1b  
**Synonyms:** IL-1F2  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >RR211023 representing NM\_031512  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCAACTGTCCCTGAACTCAACTGTGAAATAGCAGCTTTCGACAGTGAGGAGAATGACCTGTTCTTTG  
AGGCTGACAGACCCCAAAGATTAAGGATTGCTTCCAAGCCCTTGACTTGGGCTGTCCAGATGAGAGCAT  
CCAGCTTCAAATCTCACAGCAGCATCTCGACAAGAGCTTCAGGAAGGCAGTGTCACTCATTGTGGCTGTG  
GAGAAGCTGTGGCAGCTACCTATGTCTTGCCCGTGGAGCTTCCAGGATGAGGACCAAGCACCTTCTTTT  
CCTTCATCTTTGAAGAAGAGCCCGTCTCTGTGACTCGTGGGATGATGACGACCTGCTAGTGTGTGATGT  
TCCATTAGACAGCTGCACTGCAGGCTTCGAGATGAACAACAAAAATGCCTCGTGTCTGTGACCCATGT  
GAGCTGAAAGCTCTCCACCTCAATGGACAGAACATAAGCCAACAAGTGGTATTCTCCATGAGCTTTGTAC  
AAGGAGAGACAAGCAACGACAAAAATCCCTGTGGCCTTGGGCTCAAGGGGAAGAATCTATACCTGTCCCTG  
TGTGATGAAAGACGGCACACCCACCCTGCAGCTGGAGAGTGTGGATCCCAACAATACCCAAAGAAGAAG  
ATGGAAAAGCGGTTTGTCTTCAACAAGATAGAAGTCAAGACCAAGTGGAGTTTGAGTCTGCACAGTTCC  
CCAAGTGTACATCAGCACCTCTAAGCAGAGCACAGACCTGTCTTCTAGGAAACAGCAATGGTCGGGA  
CATAGTTGACTTCACCATGGAACCCGTGTCTTCC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RR211023 representing NM\_031512  
Red=Cloning site Green=Tags(s)

MATVPELNCEIAAFDSEENDLFFEADRPQKIKDCFQALDLGCPDESIQLQISQQHLDKSFRKAVSLIVAV  
 EKLWQLPMSCPWSFQDEDPSTFFSFI FEEEPVLCDSWDDDDLLVCDVPIRQLHCRLRDEQQKCLVLSDDPC  
 ELKALHLNGQNISQQVVFMSFVQGETSNDKIPVALGLKGKKNLYLSCVMKDGTPTLQLESVDPKQYPKKK  
 MEKRFVFNKIEVKTKVEFESAQFPNWIISTSAEHRPVFLGNSNGRDIIVDFTMPEVSS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/ja3076\\_f06.zip](https://cdn.origene.com/chromatograms/ja3076_f06.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF

**ACCN:** NM\_031512

**ORF Size:** 804 bp


**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
- RefSeq:** [NM\\_031512.2](#), [NP\\_113700.2](#)
- RefSeq Size:** 1339 bp
- RefSeq ORF:** 807 bp
- Locus ID:** 24494
- Cytogenetics:** 3q36
- MW:** 30.7 kDa
- Gene Summary:** an inflammatory cytokine [RGD, Feb 2006]

**Product images:**



Circular map for RR211023