

Product datasheet for **RR211030**

Sgca (NM_001107039) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Sgca (NM_001107039) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Sgca
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR211030 representing NM_001107039
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTAGCAGCCGCAACTTGGATTCTCTCCTTGCAAGTCTCCTGGCAGGCTGAGGGACAGCAAGGCC
AGCAGACAACTTTACACCACTTGTGGTCTGTGTTTGTGCATCCTTTGGAACATGCTACCTTCTGCG
CCTTCCAGAACACATTGCGGTGCCACCCACTGTCCGACTCACCTACCATGCTCACCTCCAGGACATCCA
GACCTGCCAGGTGGCTGCACTACACACAACGCAGTCCCTATAACCCCTGGCTTCTCTATGGCACTCCCA
CTCCGGAAGATCGTGGGCACCAAGTCATTGAGGTCACAGCCTACAATCGTGACAGCTTTGACACCACTAG
ACAGAGGCTGCTGCTGCTGATTGGGGACCCGAAGGTCCAGGTTACCATAACCAAGCTGAATTCCTGGTG
CGTAGCCATGATGTGGAGGAGGTCCTGCCCTCCACACCTGCCAACCCTTCTCACCGCCTTGGGGGAC
TCTGGGAGCCAGGGGAGCTCCAGTTGCTCAACATCACCTCTGCCTTGACCGAGGAGGCCGAGTCCCTCT
TCCTATTGAGGGCCGGAAGGAAGGGTTTACATTAAGGTGGGCTCTGCCACACCTTTTCCACCTGCCTG
AAGATGGTGGCATCACCTGACAGCTATGCCGTTGTGCCAGGGACGCCTCCACTATTGCTCTGCTATG
ACACCTTGGCACCTCACTCCGTGTCGACTGGTGAATGTATCTCTGGTGGACAAGTCAGTGCCAGAGCC
CCTGGATGAGGTGCCTACTCCTGGCGATGGATCTTGGAGCACGACCCATTCTTCTGCCCGCCACTGAG
GCCACGGGCCGAGACTTCTTAACAGATGCCTTGGTGACCTCTTAGTGCCCGTGCTGGTGGCTCTGCTGC
TCAACCTGTTGCTGGTTTACATTATGTGCTGTCCGCGTGAAGGACGGCTGAAGAGAGACATGGCCACCTC
TGACATCCAGATGTTTACCACCTGCTCCATCCACGGGAATACGGAAGAGCTTCGGCAGATGGCAGCCAGC
CGCGAGGTGCCCCGGCTCTATCCACCTTGCCCATGTTCAACGTGCGTACAGGAGAGCGGTTACCTCCCC
GAGTAGACAGCGCACAGATGCCCTCATCTGGACCAGCAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

ORF Size:	1161 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001107039.1 , NP_001100509.1
RefSeq Size:	1693 bp
RefSeq ORF:	1164 bp
Locus ID:	303468
Cytogenetics:	10q26
MW:	43.2 kDa