

Product datasheet for **RR211040L3V**

Angptl4 (NM_199115) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Angptl4 (NM_199115) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Angptl4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_199115
ORF Size:	1215 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211040).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_199115.2 , NP_954546.1
RefSeq Size:	1871 bp
RefSeq ORF:	1218 bp
Locus ID:	362850
UniProt ID:	Q6TMA8
Cytogenetics:	7q13
Gene Summary:	a circulating protein which causes an increase in plasma very low density lipoprotein by inhibition of lipoprotein lipase activity [RGD, Feb 2006]



[View online »](#)