

## Product datasheet for RR211040L4

### Angptl4 (NM\_199115) Rat Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Angptl4 (NM_199115) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Angptl4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211040).
Restriction Sites:	Sgfl-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.



[View online »](#)



RefSeq ORF: 1218 bp

Locus ID: 362850

UniProt ID: [Q6TMA8](#)

Cytogenetics: 7q13

Gene Summary: a circulating protein which causes an increase in plasma very low density lipoprotein by inhibition of lipoprotein lipase activity [RGD, Feb 2006]