

Product datasheet for **RR211113**

Gpr84 (NM_001109509) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gpr84 (NM_001109509) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Gpr84
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR211113 representing NM_001109509 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTGGAACAGCTCAGATGACAACCTTCCTGCTACCATGAGTCTGTATTGGGCTATCGATACTTTGCAG
TTATCTGGGGCATGGTAGTGGCTGCAACAGGCACCGTGGGCAATGTGCTCACCCCTGTTGGCCTTGCCAT
CCGTCCCAAACCTCGAACCCGTTTCAACCTGCTCATTGCCAACCTCACCCCTGGCTGATCTACTCTACTGC
ACGCTCCTGCAGCCTTTCTCCGTGGACACATACCTCCACCTCCATTGGCGCACCGGCCCATCTTCTGTA
GAATATTCGGACTCCTCCTCTTTACTTCCAATTCTGTCTCCATTCTTACCCTCTGTCTCATTGCTCTAGG
ACGCTACCTTCTCATTGCCACCCCTAAGCTCTTCCCCAGGTTTTTCAGTGCCAAGGGGATCGTGCTGGCA
CTAGTGGGAGCTGGGTTGTGGGGTGACCAGCTTTGCCCCCTCTGGAATGTTTATGTCTTGGTGCCAG
TTGTCTGCACCTGCAGCTTTGACCGGTGCGAGGCCGGCCTTACACCACCATCCTCATGGGCATCTTCTT
TGTGGTTGGGCTCAGCAGCGTGGCGTCTTCTACTGCCTCATCCACCGCCAAGTGAAGCGTGGGCTCGA
GCGCTGGACAAAATATGGGCTGCAGGAGGCCAGCATGCGCTCCATCAGGTGTCTGGACACATGAAGCTG
TGCCAGGCCACTTCCAGGAGCTAGACAGCGGGCTTGCATCAAGAGGTCCCAGCGAAGGGATTTTCATCTGA
GCCAGTCAGTGCTGCGACGACACAGACCCCTGGAAGGTGATTCGTCAGAAGCGGGGACCAGGGCATGAGA
AAGGCAGCTCAGCAGATCTCAGAGAGAAGCCTTCCAGAAGTGCATCGCAAGACTGGAGGAGCTGCAGGAG
CACGCAGAGCCACGGATGCACCATCGGAGTTCGGGAAGGTGACCCGATATGTGCTTTCAGTGTTCCTTTG
CTTCGCTCCTCAGCTACATCCCTTTCCTGCTGCTCAACATCTGGACGCCAGGGGCCGCTCCACAGGTA
GTGCATATGGTTGCTGCCAACCTCACCTGGCTCAACAGCTGCATCAACCTGTGCTCTATGCAGCCATGA
ACCGCCAGTTTCGCCAGGCTTATGGCTCCATCCTGAAACGCGGGCCACAGAGTTTCCGACGGTCCAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

ORF Size:	1188 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001109509.1 , NP_001102979.1
RefSeq Size:	1493 bp
RefSeq ORF:	1191 bp
Locus ID:	688730
Cytogenetics:	7q36
MW:	43.7 kDa