

## Product datasheet for **RR211163L4V**

### Egflam (NM\_001108938) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Egflam (NM_001108938) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Egflam
Synonyms:	RGD1306592
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001108938
ORF Size:	3015 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211163).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001108938.2</a> , <a href="#">NP_001102408.2</a>
RefSeq Size:	3714 bp
RefSeq ORF:	3018 bp
Locus ID:	365691
UniProt ID:	<a href="#">B4F785</a>
Cytogenetics:	2q16



[View online »](#)

**Gene Summary:**

Involved in both the normal retinal photoreceptor ribbon synapse formation and physiological functions of visual perception. Necessary for proper bipolar dendritic tip apposition to the photoreceptor ribbon synapse. Promotes matrix assembly and cell adhesiveness (By similarity).[UniProtKB/Swiss-Prot Function]