

Product datasheet for RR211168

Tmem5 (NM_001024759) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tmem5 (NM_001024759) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Tmem5
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR211168 representing NM_001024759
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATCGCGCTGACGCGGACCCGGCTCTGCTCGCTCCTCGTCGCCCTGTACTGCCTCTTCTCCATCTACGCCG
CCTACCACGTCTTCTCGGGCGCCGCCGGCCGGCGCGACCTCTCGGGCTCCAGGAAGACAGC
GGCGCAAGCGAAGGAGAGCGCGCGCGGTAAGCAGTCCGCGTTGGAGAGTGAAGAATGGAATCCTTGG
GAAGGAGATGAAAAAACGAGCAGCAACACAGAATTTAAACTGGCCTGCAAATCTTAAACAAATCCACAA
AAGAAAAATAGAGCATCGTGTGCAAATCTGGGCAAAAGCTGCCATCGGTCTCTATCTCTGGAAACATGT
TTTTGAAGGCACTCTCGCTCCTGCTGATGTGACTGCTCAGTGGAGGGAAGGACAGTCAAGTGTAGGAAGA
ACACATTACAGCTTCATCACTGGCCAGCTGTAGTTCCCGGTACTTCTCCATAGATGTGGACAATGTGG
TCCTTGTTTTAAATGGAAGAGAAAAAGCAAAGATCTTCCATGCCACCCAGTGGTTAACTTATGCCCAGAA
TTAATAACAAACCCAAAACTGCAGCATCTGGCTGTTGTCTTGCTTGGAGATGAGCACTGTGATAATGAC
TGGATAATGCAGTTCCTCAGAAGAAATGGAGGCTTCGTGGACCTACTTTTCATAACATATGACAGCCCTT
GGATTAATGGTGCAGATGTCCTCAGTGGCCGTTAGGTGTAGCAACATACAGGAAGTTTCTGTGGTTGA
AGCCAGCTGGTCCATGCTGCATGATGAGAGACCTATATATGTAATTTCTTAGGGACTATCTATGAAAT
TCATCAAGACAAGCACTAATGAACATTTTAAAACAAGATGGAAATGATAAACTTTGTTGGTTTTCTGCAA
GAGAACAGTGGCAGCCTCAGGAGACAAATGAAAGCCTTAAAAATACCAAGATGCTTTGCTTCATAGTGA
CCTCACATTGTGCCCTGTAGGAGTGAACACTGAATGTTACAGGATCTATGAGGCCGTCTCTTTGGCTCC
ATTCCTGTGGTAGAAGACGTGATGACAGCTGGCCACTGTGAAACACATCCAGTCAGCACAGCGCTCCTC
TGCAGTACTCAAGGCCATGGGGCCCCCTTCACTTCAAGAAGTGGAAAGGAGCTCCCTGCTATTTT
AGAAAAGGAAAAAATAATCTTACAAGAAAAGATTCAAAGAAGAAAAATGTTACTTTCATTGGTACCAA
CATTTCAAACAGAACTAAAGTGAAATTTACTAAAATTTAGAAAGTTCATTTTTATGAATAAAG
GT

ACGCGTACGCGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR211168 representing NM_001024759
 Red=Cloning site Green=Tags(s)


MRLTRTRLCSELLVALYCLFSIYAAYHVFFGRRRRPPGATSRGSRKTAQAQAKERRGRGKQSALESEEWNPW
 EGDEKNEQQHRIKTGLQILNKSTKEKIEHRVQIWGKAAIGLYLWEHVFEGLAPADVTAQWREGQSVVGR
 THYSFITGPAVVPGYFSIDVDNVVLVLNGREKAKIFHATQWLTYAQNLIQTQKLOHLAVVLLGDEHCDND
 WIMQFLRRNGGFVDLLFITYDSPWINGADVLQWPLGVATYRKFVVEASWSMLHDERPYICNFLTGIYEN
 SSRQALMNILKQDGNKLCWVSAREQWQPQETNESLKNYQDALLHSDLTLCPLVGNTECYRIYEACSFGS
 IPVVEDVMTAGHCGNTSSQHSAPLQLLKAMGAPFIFIKNWKELPAILEKEKTIILQEKIQRKMLLHWYQ
 HFKTELKWKFTKILESSFFMNNKG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_001024759
ORF Size:	1332 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001024759.1 , NP_001019930.1
RefSeq Size:	1549 bp
RefSeq ORF:	1335 bp
Locus ID:	299841
Cytogenetics:	7q22
MW:	51.2 kDa