

## Product datasheet for **RR211185L4V**

### **Pigl (NM\_138901) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Pigl (NM_138901) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pigl
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_138901
ORF Size:	756 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211185).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_138901.2</a> , <a href="#">NP_620256.1</a>
RefSeq Size:	1410 bp
RefSeq ORF:	759 bp
Locus ID:	192263
UniProt ID:	<a href="#">O35790</a>
Cytogenetics:	10q23
Gene Summary:	catalyzes the N-deacetylation of N-acetyl-D-glucosaminylphosphatidylinositol in the second step of GPI anchor biosynthesis [RGD, Feb 2006]



[View online »](#)