

Product datasheet for RR211195L3

Dgka (NM_080787) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Dgka (NM_080787) Rat Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Dgka
Synonyms:	Dagk1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211195).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.



[View online »](#)

Locus ID: 140866
UniProt ID: [P51556](#)
Cytogenetics: 7q11
Gene Summary: may mediate myelin formation and metabolism; expressed specifically in oligodendrocytes [RGD, Feb 2006]