

Product datasheet for **RR211224**

Rhot1 (NM_001107026) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rhot1 (NM_001107026) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Rhot1
Synonyms:	MGC156790
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RR211224 representing NM_001107026
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCGTGGCGGCCGAGTGGCCCTCGGAGCCTCCGACATGAAGAAAGACGTGCGGATCCTGCTGGTGG
 GAGAACCTAGAGTTGGGAAGACATCACTGATTATGTCTCTGGTCAGTGAAGAATTTCCGGAAGAGGTTCC
 TCCTCGGGCAGAAGAAATCACCATTCCAGCTGATGTCACCCCTGAAAGAGTTCCAACACACATTGTTGAT
 TATTCAGAGGCAGAACAGAGTGATGAACAACCTTCATCAAGAAATATCACAGGCTAATGTCATCTGCATAG
 TTTATGCTGTTAAACAACAAACATTCTATTGATAAGGTAACAAGCCGATGGATTCTCTCATAAATGAAAG
 AACAGACAAAGACAGCAGGCTGCCTTTGATATTGGTTGGGAACAAATCTGACCTGGTGGAGTATAGTAGT
 ATGGAGACCATCCTTCCAATTATGAACCAATATACAGAAATAGAAACCTGTGTGGAGTGTTCAGCAAAAA
 ACCTAAAGAATATATCAGAACTGTTCTATTATGCACAGAAAGCTGTTCTTCATCTACAGGGCCTTTGTA
 CTGCCAGAGGAGAAAGAGATGAAACCAGCTTGATAAAAAGCCCTCACTCGTATATTTAAAAATATCTGAC
 CAAGATAATGATGGCACTCTGAATGATGCTGAGCTCAACTTCTTTAGAGAAATTTGTTTCAACACTCCTT
 TAGCTCCCAAGCTCTGGAAGATGTCAAGAATGTGGTTAGAAAGCACTTAAGTGACGGTGTGGCTGACAG
 TGGGCTGACACTCAGAGGTTTTCTTTTTACATACACTTTTTATCCAGAGGGGGAGGCATGAGACCACT
 TGGACTGTGCTTCGACGATTTGGTTATGATGATGATCTGGATCTGACGCCGAGTACTTATTTCCCTGC
 TAAAAATACCTCCTGATTGCACTACTGAACTAAATCATCATGCATATTTGTTTCTCAAAGCACCTTTGA
 CAAACATGATCTGGACAGAGACTGTGCTTTGTACCTGATGAACTCAAAGATTTATTTCAAGTTTTCCCT
 TACATACCTTGGGACCAGATGTAATAACACAGTTTGTACAAATGAAAGCGGCTGGATAACCTACCAGG
 GATTCTTTCCAGTGGACGCTCACAACCTACTTAGATGTACAGGTTGTCTGGAGTATCTGGGCTATCT
 AGGCTATTCATACTGACGGAACAAGAGTCTCAAGCTTCAGCCATCACAGTAACAAGAGATAAAAAGATT
 GACCTGCAGAAAAACAGACTCAGAGAAATGTGTTGAGATGTAATGTGATTGGAGTGAAGGCTGTGGAA
 AAAGTGGAGTTCTTCACTCTCTCTTGGAAAGAACTAATGAGGCAGAAGAAAATTCGTGATGATCATAA
 ATCCTATTATGCAATTAACACTGTCTATGTATACGGACAAGAGAAATACTTGTGCTACATGATATCTCA
 GAGTCGGAGTTTCTTACTGAGGCTGAGATCGTTTGTGATGTTGTATGCCTGGTATATGATGTCACCAATC
 CCAAGTCATTTGAATACTGCGCCAGGATTTTTAAGCAACATTTTATGGACAGCAGAATTCCTGCTTAAT
 TGTAGCTGCAAAGTCAGACCTGCATGAAGTCAAGCAAGAACACAGTGTCTCACCTACTGATTTCTGCAGG
 AAACACAAAATGCCTCCGCTCAAGCCTCACTTGTAACTGCTGATGCGCCAGTAAGGATATCTTTG
 TTAAGCTGACAACGATGGCCATGTACCCGACGTGACACAAGCTGACCTCAAGAGCTCAACATTCTGGCT
 TCGAGCAAGCTTTGGTGTACTGTCTTTCAGTGTGGGCTTTGCTATGTACAGGCACTATTGAAACAG
 CGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR211224 representing NM_001107026
 Red=Cloning site Green=Tags(s)

```
MRAGRVRPLRASDMKKDVRILLVGEPRVGKTSLSIMSLVSEEFPEEVPPRAEEITIPADVTPERVPTIIVD
YSEAEQSDQLHQEISQANVICIYAVNNKHSIDKVTSRWIPLINERTDKDSRLPLILVGNKSDLVEYSS
METILPIMNQYTEIETCVECSAKNLKNISELFYYAQKAVLHPTGPLYCPEEKEMKPACIKALTRIFKISD
QDNDGTLNDAELNFFQRICFNTPLAPQALEDVKNVVRKHLSDGVADSGLTLRGFLFLHTLFIQRGRHETT
WTVLRRFGYDDDLDTPEYLFPLLKIPPDCTTELNHHAYLFLQSTFDKHDLDRCALSPDELKDLFQVFP
YIPWGPVDVNTVCTNESGWITYQGFLSQWTLTTYLDVQRCLYLYGLGYSILTEQESQASAITVTRDKKI
DLQKKQTQRNVFRCNVIGVKGCGKTGVLQSLGRNLMRQKKIRDDHKSYYAINTVYVYGQEKYLLLDHDS
ESEFLTEAEIVCDVVCLVYDVTNPKSFEYCARIKQHFMDSRIPCLIVAASDLHEVKQEHSVSPTDFCR
KHKMPPPQAFTCNTADAPSKDIFVKLTMMAMYPHVTQADLKSSTFWLRASFGATVFAVVGAFMYRALLKQ
R
```


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_001107026

ORF Size: 1893 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001107026.2](#), [NP_001100496.1](#)

RefSeq Size: 2941 bp
RefSeq ORF: 1896 bp
Locus ID: 303351
Cytogenetics: 10q25-q26
MW: 72.2 kDa