

## Product datasheet for **RR211250L4V**

### TPP1 (NM\_031357) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TPP1 (NM_031357) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	TPP1
Synonyms:	Cln2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031357
ORF Size:	1689 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211250).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_031357.1</a> , <a href="#">NP_112647.1</a>
RefSeq Size:	2485 bp
RefSeq ORF:	1692 bp
Locus ID:	83534
UniProt ID:	<a href="#">Q9EQV6</a>
Cytogenetics:	1q32
Gene Summary:	human homolog catalyzes the removal of N-terminal tripeptides from polypeptides, also has a minor endoprotease activity [RGD, Feb 2006]



[View online »](#)